

Product Data Sheet

Product Name: Deluxe CPR Mask

Item No. : 1800061

Product Character:

The shield offers you personal safety and it prevents physical contact with victim face during rescuing breathing. Anti-re flux vale opens with a butterfly design to prevent inverse flow from victims breath or body fluids through the airway.

Product Pictures:



Specification :

Item	Details
Product Name	Deluxe CPR Mask
Item No.	1800061
Place of Origin	China
Product Classification	Class I
Material	EVA
Total Pieces	1pc
Mask Size	20*28cm
Color	Transparent
Brand	Ever Ready First Aid
Application	For cardiopulmonary resuscitation, protective barrier should be used
Type	disposable face shields
Certificate	CE, ISO,FDA
Packing	1pcs/PE Bag
Export Carton Details	60*42*36 cm 1000pcs/ctn G.W.: 9.8KG
H.S.Code	3926909090

info@everreadyfirstaidusa.com

www.everreadyfirstaidusa.com



Instruction for use:

Cardiopulmonary resuscitation (CPR)

Warning:

- 1、 This product is used by people who have been trained by CPR.
- 2、 Once the patient begins to breathe spontaneously, the rescuer must remove the product from the patient.
- 3、 The product is disposable, respiratory valve blockage or potential blockage caused by stomach or saliva, no further use. Do not pasteurize, boil, steam or gas autoclave this device.