



MEDITAC BASIC PREMIUM BLEEDING CONTROL PACK

Feat. C.A.T. Tourniquet, Israeli Bandage and NAR Compressed Gauze Dressing

PRODUCT INCLUDES:

- Basic premium bleeding control pack is a single-use kit for traumatic blood loss.
- Non-resealable package with Quick Tear feature for easy access.
- C.A.T. Tourniquet has multiple uses for injuries, minor to life-threatening.
- Israeli emergency bandage for hemorrhage control with pressure applicator.
- NAR Compressed Gauze Dressing includes unique crinkle weave for fluid absorption.
- Mylar Emergency Blanket is durable and double-sided to preserve body heat in emergency or extreme conditions.
- Nitrile Gloves-size L, tear-resistant.
- 7.5" Shears are sturdy for cutting tough cloths and bandages.
- Permanent Marker which marks and stands out on almost every surface.



Item Dimensions: 6" x 2.5" x 9"

SKU: MTBCBPB

844-300-3028

info@meditackits.com

300 Liberty Ave, Brooklyn, NY 11207

Meditac

@meditackits

@meditackits