

Dixie Ems Deluxe Aneroid Sphygmomanometer Blood Pressure Set W/Adult Cuff, Carrying Case and Calibration Tool

Reorder	Color	Count
143401	Navy	25/case
143401BK	Black	25/case
143401PK	Pink	25/case
143401PUR	Purple	25/case

Dependable instrument for EMS, hospital or home use. It features a navy nylon cuff with range markings to ensure proper cuff size selection and an artery label to facilitate correct positioning.

LATEX FREE BLOOD PRESSURE CUFF

It can be operated easily with a single finger or thumb. Made of a durable nylon that secures comfortably with a high quality Velcro strip.

HIGH CONTRAST & EASY TO READ DIAL

Allows the user to read blood pressure measurements easily and accurately.

ADULT SIZED

Nylon cuff features a size range from 25.4cm to 40.6cm, ideally sized for adult arms between 10" and 16" in circumference.

DELUXE END VALVE

Designed to reduce dust build up.

ZIPPERED CARRYING CASE

Easily store and transport when not in use with a secured zippered carrying case.

